


The 91439A Flexport® interface integrates data provided from Pulsion PiCCO₂ monitor into Spacelabs Healthcare XPREZZON®, Qube®, or Ultraview SL™ monitors. The Flexport interface provides bedside and central station display of third party device data and alarms, and provides tabular and graphical trending of cardiac data. The Flexport interface provides Remote View and Alarm Watch capabilities from the central and remote monitors. The data provided from the Flexport interface is made available through the ICS interface to clinical information systems. All compatible XPREZZON, Qube, and Ultraview SL monitors support 91439A Flexport interface.


Continuous View

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Split Screen View

Numeric Parameter key

Features

Parameters

The 91439A Flexport interface counts as one parameter when computing parameter capacity for monitors.

Product Specifications

Physical Dimensions

Height	2.5 cm (1 in)
Width	13.5 cm (5.3 in)
Depth	4.6 cm (1.8 in)
Weight	0.11 kg (4 oz)

Environmental Requirements

Operating

Temperature	10° to 40° C (50° to 104° F)
Humidity	95% (non-condensing)
Altitude	0 to 2,000 m (0 to 6,562 ft)

Storage

Temperature	-40° to 75° C (5° to 167° F)
Humidity	95% (non-condensing)
Altitude	0 to 12,192 m (0 to 40,000 ft)

Monitor Compatibility

Monitors Supported

XPREZZON®	XPREZZON (91393)
Qube®	Qube (91390)
Ultraview SL™	SL2200 (91367), SL2400 (91369), SL2600 (91370), SL2700 (91387-27), SL2800 (91387-28), SL3800 (central 91387-38), SL3900 (central 91387-39)

Monitors Not Supported

UCW® Classic (90385), Ultraview 1700 Classic (90387), Ultracare SLP100, PC2 (90305), PC1 (90303), PCX (90308), PC Scout® (90309), mCare (91220), Maternal Obstetrical Monitor (MOM™) (94000)

Note:

Spacelabs Healthcare confirms that the Flexport software versions shown in the table that follows are verified for use with the software or protocol versions listed for the source devices.

Company, Product, or Source Device	PiCCO ₂ Device Version	Spacelabs Healthcare Flexport Software Version	Flexport Model Number
PiCCO ₂	3.1.0.8A	1.00.00-013	91439A

Note:

The trending function is unavailable for the XPREZZON and qube monitors.

Flexport Interface

Note:

Each interfaced device must have an appropriate RS-232 interface and operating software of a specific level to function with Flexport interfaces. Contact the device manufacturer to confirm compatibility before ordering.

SDLC Cable	
Qube	P/N 012-0152-xx
XPREZZON	P/N 012-0555-xx
SL2200/SL2400/SL2600	P/N 012-0152-xx
SL2700/SL2800	P/N 012-0555-xx
SDLC Terminator	P/N 012-0507-02
Transition Connector	P/N 131-2391-xx

Data Shown

Label	Name
Alarm Status	Bell shows when alarms are turned on.
CI	Cardiac Index
CO	Cardiac Output
CFI	Cardiac Function Index
CPI	Cardiac Power Index
ELWI	Estravascular Lung Water Index
GEDI	Global End-Diastolic Volume Index
GEF	Global Ejection Fraction
PiCCO	PiCCO2 Flexport key

Label	Name
PCCI	Pulse Contour Cardiac Index
PCCO	Pulse Contour Cardiac Output
PVPI	Pulmonary Vascular Permeability Index
ScvO2	Central Venous Oxygen Saturation
SVI	Stroke Volume Index
SVRI	Systemic Vascular Resistance Index
SVV	Stroke Volume Variation
TD Values	Thermo Dilution Values Time and Date
VO2I	Oxygen Consumption Index

Accessories

Connectors	
Female	6-pin modular jack
Male	9-pin SDLC connector
Cable Length	
SDLC	30.48 cm (12 in)
RS-232	243.8 cm (96 in)

Maximum SDLC Cable Length
The total length of all SDLC cables from the power source (monitor or module housing) to the Flexport system interface should not exceed 1.8 m (6 ft).
When linking other Spacelabs Healthcare devices, the Flexport system interface must be the closest device to the monitor or module housing.

Documentation

CD-ROM Part Number *91439A Flexport Interface, PiCCO₂ Operations and Service CD-ROM*
(P/N 084-1601-xx)

Supplies and Accessories *Spacelabs Healthcare Supplies and Accessories Catalog*
(sa.spacelabshealthcare.com)

Regulatory Approvals

CE marked in accordance with the Medical Device Directive 93/42/EEC.

Please refer to <https://www.spacelabshealthcare.com/about-us/patents-trademarks> for a full listing of Spacelabs Healthcare trademarks. Other brands and product names used herein are trademarks of their respective owners.

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