da Vinci® Surgical System User Manual

Table of Contents

Intuitive Surgical, Inc. v

P/N 550120-02 Rev XB

How to Use This Manual	xi
General Information	xii
Contact Information	xii
General Information for U.S. Only	xiii
Essential Prescribing Information	xiii
Indications for Use Statement	xiii
Professional Instructions for Use	xiii
General Precautions, Warnings, and Contra indications	xv
Conversion to Non-Minimally Invasive Technique	xv
Endoscopic Procedure Precautions	xv
High Frequency Electrosurgery Precautions	xvi
Installation and Service Precautions	xvii
Transportation and Storage Precautions	xvii
International Electrotechnical Commission (IEC) 601 General Information	xviii

Instrument and Endoscope Isolation	xviii
Accessory Equipment Interconnection	xviii
General Requirements	xix
Power Requirements	xix
Environmental Conditions: Operating	xix
Environmental Conditions: Storage and Transport	xix
Electromagnetic Compatibility	xix
Battery Disposal Information	xx
Chapter 1: Introduction to the da Vinci System	1-1
Masters	1-2
3D Display System and Stereo Viewer	1-2
Controls and Indicators	1-2
Surgical Cart Overview	1-4
The InSite® Vision System Overview	1-6
Chapter 2: Setting up the da Vinci System	2-1
Setting Up the Operating Room	2-1
Surgeon Console Positioning	2-1
Main Circuit Breakers	2-2
Surgical Cart Positioning	2-2
Vision Cart Positioning	

2-3
 Cable Connection Guide 2-3
 Start-Up
 2-5
 Surgeon Console and Surgical Cart
 2-5

Table of Contents

vi Intuitive Surgical, Inc.

P/N 550120-02 Rev XB

Chapter 3: Draping the da Vinci System 3-1

Prepare the Sterile Accessories
 3-1

Sequence of Draping
 3-2

Guidelines for Instrument Arm Draping
 3-2

Attaching Sterile Adapters to the Drapes 3-2

Attaching Cannulae to the Instrument Arms 3-4

Camera Arm Draping Instructions
 3-5

Draping Instructions For the Stereo 3D Camera
 3-7

The InSite® Vision System Set-Up
 3-8

Camera and Endoscope Calibration
 3-9

Operating Modes
 3-10

Stand By Mode
 3-10

Ready Mode
 3-10

Follow Mode
 3-11

Emergency Stop
 3-12

Instrument Arm Sterile Adapters Compatibility	3-13
Chapter 4: Use During a Procedure	4-1
Selecting Port Locations	4-1
Using the Surgeon's Console	4-2
Adjusting the Console for Comfort	4-2
Working with the Video Image	4-2
Selecting the Scope Angle	4-2
Switching Display from 3D to 2D	4-3
Focusing the Video Image	4-3
Repositioning the Camera/Scope Image	4-3
Camera Orientation Line	4-4
Using the Auxiliary Instrument Arm	4-4
Operating the Masters	4-5
Repositioning the Masters	4-5
Scaling	4-5
Working with the Surgical Cart	4-6
Cannula Insertion and Connection	4-6
Selecting the Cannula Length	4-6
Camera Cannula Connection to Camera Arm	4-7
Instrument Arm Cannula Insertion and Connection	4-7
Option 1: Instrument Arm Cannula Not Attached To The Instrument Arm	4-7
Option 2: Instrument Arm Cannula Attached To The Instrument Arm	4-8

Instrument Arm and Camera Arm Clutching 4-8

Working with Instruments 4-9

Instrument Placement and Removal 4-9

Manually Inserting Instruments 4-10

Guided Tool Change 4-10

Removing Instruments 4-11

Using the Blade Protector 4-12

Table of Contents

Intuitive Surgical, Inc. vii

P/N 550120-02 Rev XB

Instrument Grip Release 4-13

Tool Reader Mode 4-14

Scalpel Blade Insertion Instructions 4-14

Instructions for Loading and Firing of Clip Applier Instruments 4-15

Loading and Firing Instructions 4-15

Fluid Leakage Precautions 4-16

The InSite® Vision System 4-17

Camera/Endoscope Attachment and Removal 4-17

Chapter 5: Shutting Down the System 5-1

Preparing for Shutdown 5-1

Instrument Usage Summary Display 5-2

Shutdown Process 5-3

Conversion to Open Surgery 5-4

Chapter 6: Understanding the da Vinci

Components	6-1
Surgeon's Console	6-1
Controls and Indicators	6-1
User Interface Panel (UIP) and User Switch Panel (USP)	6-2
User Switch Panel – Mode Control	6-3
SYSTEM (On/Off)	6-3
EMERGENCY STOP	6-3
STAND BY	6-4
READY	6-4
User Switch Panel – Status Indicators	6-4
User Interface Panel – Controls and Indicators	6-5
System Status	6-5
SCALING	6-5
SCOPE ANGLE	6-5
DISPLAY	6-6
CANNULA	6-6
CAMERA CALIBRATION	6-7
SELECT MASTER	6-7
Disabling an Auxiliary Instrument Arm	6-8
Viewing Auxiliary Arm Information	6-8
Other Controls on the Surgeon Console	6-10
Head Sensor	6-10
Viewer Height	6-10
Intra-ocular Spacing Adjustment	6-11

Audio Intercom System Volume Adjustment 6-11

Brightness and Contrast Adjustment 6-11

Emergency Power Off (EPO) 6-11

Icons 6-12

Foot Switches 6-14

Clutch 6-14

Table of Contents

viii Intuitive Surgical, Inc.

P/N 550120-02 Rev XB

Camera 6-14

Focus 6-15

Unlabeled 6-15

Coag 6-15

Surgical Cart 6-15

Set-up Joints 6-15

Instrument Arms 6-16

Camera Arm 6-16

Auxiliary Arm 6-16

Moving the Auxiliary Arm (Z-axis) 6-17

Collision Detection 6-18

Excessive Force Detection 6-19

Drapes and Sterile Adapters 6-19

Instruments 6-19

Instrument Identification and Description

6-19	
Limited Reuse and Expiration Indication	6-22
Chapter 7: Cleaning and Maintenance	7-1
The da Vinci System	7-1
Surgeon Console	7-1
Surgical Cart	7-1
Instruments	7-1
Accessories	7-1
All Other Maintenance	7-3
Chapter 8: System Intelligence	8-1
Fault and Emergency Stop State Philosophy	8-1
Non-recoverable Faults	8-2
Recoverable Faults	8-2
Loss of Power/Reserve Power	8-3
Troubleshooting	8-3
Audible Alarms	8-3
Viewing the Error Log	8-4
Symbols and Icons	8-4
Chapter 9: Glossary of Terms	9-1
Appendix A: Clinical Appendix	A-1

Clinical Study A-1

Study Endpoints A-1

Study Design: Enrollment, Assessment and Randomization A-1

Randomization Results A-2

Complications/Adverse Events A-3

Procedures Converted A-3

Table of Contents

Intuitive Surgical, Inc. ix

P/N 550120-02 Rev XB

Outcomes Analysis A-4

Device Failures and Replacements A-6

Risk/Benefit Analysis A-6

Safety A-6

Effectiveness A-6

Table of Contents

x Intuitive Surgical, Inc.

P/N 550120-02 Rev XB

Intuitive Surgical, Inc. xi

P/N 550120-02 Rev XB

How to Use This Manual

The da Vinci® Surgical System User Manual and the InSite® Vision System User Manual are intended to be used together to support the da Vinci System. The da Vinci® Surgical System User Manual is organized to provide you with the steps you need to get things done first, then provide reference information should you decide to learn more about each task or component. Chapters 2-5 provide the basic instructions for setting up, using and shutting down the system. Chapters 6-8 contain detailed reference information about the components, how to clean and maintain them, and information on how the system intelligence functions. The final sections

contain a glossary, a clinical appendix describing results from clinical trials, and an index. A brief description of each chapter is provided below.

Front Matter (the chapter you are reading now)

This section contains regulatory and safety information that must be read by every user of the da Vinci System. Examples include details on general precautions, precautions specific to procedures, power information, and so on.

Chapter 1, Introduction to the da Vinci System

Chapter 1 provides a basic overview of the system, including high-level information on the components and a brief description of how the system is used, as well as features and benefits.

Chapter 2, Setting Up the da Vinci System for Surgery

This chapter offers step-by-step instructions for preparing the da Vinci System for surgery. The task-oriented steps include instructions for orienting and connecting the components with cabling, and how to power up the system.

Chapter 3, Draping the da Vinci System

Chapter 3 explains how to drape the system, how to connect the camera/endoscope assembly, and how to calibrate the camera and endoscope.

Chapter 4, Using the da Vinci System During a Procedure

Chapter 4 describes step-by-step use of the da Vinci System during a surgical procedure. Tasks include setting up the Surgeon's Console for the comfort of the surgeon; details on using the Surgical Cart; and instructions on how to use specific instruments. Chapter 4 also contains a brief discussion on the use of the InSite Vision System.

Chapter 5, Shutting Down the da Vinci System

This chapter explains how to prepare for a shutdown, and explains how to track instrument usage.

Chapter 6, da Vinci Component Reference

Chapter 6 describes each system component in detail. Use this chapter to get more information if the basic steps in Chapters 2-5 do not provide enough detail.

Chapter 7, Cleaning and Maintenance

The instructions for recommended care of all parts of the da Vinci System, including the instruments, are explained in detail. Suggested methods for sterilization are included as well.

Chapter 8, System Intelligence

This reference-oriented chapter contains descriptions of Faults and the Emergency Stop state, when faults can occur, and how to manage them. A listing of all system icons is provided as well.

Chapter 9, Glossary of Terms

Chapter 9 provides an explanation of terms and acronyms used in this manual.

Appendix A, Clinical Appendix

The Clinical Appendix provides specific details of clinical trials using the da Vinci System.

General Information

The Intuitive Surgical da Vinci Endoscopic Instrument Control System, Model IS1200, is in conformance with the Medical Device Directive, 93/42/EEC.

Contact Information

Manufacturer: Intuitive Surgical, Inc.
950 Kifer Road
Sunnyvale, CA 94086-5206
USA
Toll free: 1-888-868-IMIS (4647)
Phone: 408-523-2100
Fax: 408-523-1390

European Representative Intuitive Surgical, S.A.R.L.
5, Place Royale
78100 Saint-Germain en Laye
France
Phone: +33-1-39-04-26-60
Fax: +33-1-39-04-26-61

The da Vinci Endoscopic Instrument Control System (Model IS1200) is manufactured in the USA.

General Information for U.S. Only

The information in this section is intended for U.S. customers. General information continues in the next section, "General Precautions, Warnings, and Contra indications" on page xv.

Essential Prescribing Information

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Indications for Use Statement

Device Name: Intuitive Surgical® da Vinci® Endoscopic Instrument Control System and Endoscopic Instruments.

Indications for Use: The Intuitive Surgical da Vinci Endoscopic Instrument Control System (hereinafter referred to as the "da Vinci System") is

intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments including rigid endoscopes, blunt and sharp endoscopic dissectors, scissors, scalpels, forceps/pick-ups, needle holders, endoscopic retractors, stabilizers, electrocautery and accessories for endoscopic manipulation of tissue, including grasping, cutting, blunt and sharp dissection, approximation, ligation, electrocautery and suturing during general laparoscopic surgical procedures, general non-cardiovascular thoracoscopic surgical procedures, and thoracoscopically-assisted cardiotomy procedures. It is intended to be used by trained physicians in an operating room environment in accordance with the representative, specific procedures set forth in the Professional Instructions for Use.

Professional Instructions for Use

General Indication for Use: The da Vinci System is intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments including rigid endoscopes, blunt and sharp endoscopic dissectors, scissors, scalpels, forceps/pick-ups, needle holders, endoscopic retractors, stabilizers, electrocautery and accessories for endoscopic manipulation of tissue, including grasping, cutting, blunt and sharp dissection, approximation, ligation, electrocautery and suturing.

Representative Uses: The da Vinci System has been successfully used in the following procedures, among others:

- General laparoscopic surgical procedures such as cholecystectomy, Nissen fundoplication, radical prostatectomy, Heller myotomy, gastric bypass, donor nephrectomy, adrenalectomy, splenectomy and bowel resection.
- General non-cardiovascular thoracoscopic surgical procedures such as internal mammary artery mobilization.
- Thoracoscopically-assisted cardiotomy procedures such as mitral valve repair and totally endoscopic atrial septal defect closure.

WARNING: The system should be used only by surgeons who have developed adequate robotic skills to perform the tasks associated with each procedure and who have received specific training provided by Intuitive Surgical, Inc., in the use of this device.

WARNING: Performance characteristics for conduct of totally endoscopic or minimally invasive coronary artery bypass surgery using the da Vinci System have not been established.

WARNING: The clinical evaluation of the da Vinci System supporting its use for mitral valve

repair was not performed totally endoscopically. Introduction and manipulation of the endoscopic instruments were controlled by the da Vinci System through port incisions (< 1 cm) while accessory technologies, such as atrial retractor and cardioplegia line, etc., were introduced through a mini-thoracotomy. Performance characteristics for conduct of totally endoscopic mitral valve repair using the da Vinci System have not been established. Therefore, the system should be used only by surgeons who have developed adequate robotic skills to perform the tasks associated with each procedure and who have received specific training provided by Intuitive Surgical, Inc. in the use of this device.

WARNING: The friable nature of pulmonary tissue enhances the risk of vascular, bronchiolar, or other injury that will be difficult to control when using this device. Published clinical experience as well as clinical studies performed to support this marketing clearance have demonstrated that even surgeons considered expert in laparoscopy/thoracoscopy have substantial learning curves of 10 to 12 cases.¹ Therefore, the system should be used only by surgeons who have developed adequate robotic skills to perform the tasks associated with each procedure and who have received specific training provided by Intuitive Surgical, Inc., in the use of this device.

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所有直接或間接相關的內視鏡手術之禁忌症候 **Any and all relative and absolute contra indications to endoscopic surgical technique applicable to the use of conventional endoscopic surgical instruments apply to the use of the da Vinci System. General, non-procedure specific, contra indications to endoscopic surgery include bleeding diathesis, morbid obesity and pregnancy.**

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¹ Falk, et al., Total endoscopic computer enhanced coronary artery bypass grafting, Eur J Cardiothorac Surg 2000; 17: 38-45.

CAUTION: Thermal hazards may exist from high temperatures. Eye hazards may exist from the high energy radiated light by the endoscopic camera and illumination system. Only personnel having adequate training and experience with the endoscopic camera and illumination system should operate such equipment. All WARNING and CAUTION messages provided with the endoscopic camera and illumination system must be followed.

CAUTION: The force feedback associated with the da Vinci System is different from feedback experienced when using conventional instruments. As with any endoscopic procedure, the surgeon should rely on visual cues to enhance force feedback.

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CAUTION: Do not use electroSurgical equipment unless properly trained in the specific procedure being undertaken. Follow all instructions, warnings, and cautions provided with the ESU.

CAUTION: The Intuitive Surgical ElectroSurgical Instrument is designed for use with the coagulation (COAG) waveform of high frequency ESUs and a maximum peak voltage of 3kV (6kV peak-to-peak). Do not use settings on the ESU that exceed 3kV peak. Do not attempt to use the Foot Switch on the Surgeon Console with ESUs that are not compatible with the da Vinci System. Consult with your Intuitive Surgical Representative regarding compatible models.

Note: Only the COAG function is wired to the Foot Switch control.

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CAUTION: Do not clean the tip with another instrument using directional movements away from the tip.

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If the da Vinci System does cause harmful interference to other devices, the user is encouraged to try to correct the interference by one of more of the following measures:

- Reorient or relocate the affected device(s)
- Increase the separation between the da Vinci System and the affected device(s)
- Connect the da Vinci System into an outlet on a circuit different from that to which the

affected device(s) is connected

- Contact your Intuitive Surgical Representative for help

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

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The da Vinci System contains a non-spillable, lead acid battery pack module that is not user serviceable. This battery pack must be disposed of in accordance with local regulations. Contact your local Intuitive Surgical Representative for disposal information.

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
Intra-ocular Spacing Adjustment

The intra-ocular spacing adjustment knob is located beneath the Stereo Viewer, and is used to optimize the fusion of the 3D operative image. It adjusts left/right spacing of the image to match the spacing of the surgeon's eyes.

Beneath Stereo Viewer on Surgeon's Console

- A. Intra-Ocular Spacing Adjustment
- B. Brightness and Contrast Knobs for Stereo Viewer
- C. Audio Intercom System Volume Knob

Figure 6-7 Location of Stereo View Adjustment Knobs

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Excessive Force Detection

The auxiliary arm is constructed so the system becomes aware of any excessive force generated during movement of the arm along the Z-axis, and stops the movement of the arm in that situation. For example, it is possible that while raising the auxiliary arm, the arm could collide with the bottom of the operating table. The system recognizes when the arm begins to generate force on an object, and automatically stops the lift motor.

A repeated tone sounds and an icon is displayed until the force is removed. The force can be removed by moving the arm in the opposite direction to clear the fault, or simply removing the obstruction.

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Instrument Descriptions	
	<p>Long Tip Forceps</p> <p>Designed for atraumatic grasping of tissue.</p>
	<p>Round Tip Scissors</p> <p>Designed to cut tissue or suture.</p>
	<p>DeBakey Forceps</p> <p>Designed for atraumatic grasping of tissue.</p>
	<p>Small Clip Applier</p> <p>Specifically designed to apply clips to ligate vessels.</p>
	<p>Scalpel/Cautery</p> <p>Designed to hold a 15-degree scalpel blade. The blades are disposable and can be inserted into and removed from the instrument.</p>
	<p>Scalpel/Cautery</p> <p>Designed to hold a Beaver blade. The blades are disposable and can be inserted into and removed from the instrument.</p>
	<p>Large Needle Driver</p> <p>Has carbide inserts and is specifically designed to grasp and drive needles.</p>

	<p>Potts Scissors Designed for delicate cutting of tissue.</p>
	<p>Round Tooth Forceps Designed for atraumatic grasping of tissue.</p>
	<p>Black Diamond Micro Forceps Has a delicate instrument tip and is designed to grasp smaller needles as well as tissue.</p>
	<p>Cadiere Forceps Designed for atraumatic grasping of tissue.</p>
	<p>Cautery with Spatula Designed for blunt dissection and the application of cautery.</p>
	<p>EndoWrist ProGrasp™ Tissue grasper with greater grasping force than the Cadiere Forceps.</p>
	<p>Cichon Tissue Forceps Designed for atraumatic grasping of tissue. The Cichon Tissue Forceps provide similar capabilities as the DeBakey Forceps, but with a much smaller instrument profile that increases the surgeon's exposure to the surgical site.</p>

	<p>EndoWrist PreCise™ Bipolar Instrument</p> <p>This is a fully articulating instrument that provides precise electrocautery and hemostasis. Designed for atraumatic grasping of tissue and the application of bipolar cautery.</p>
	<p>Ultrasonic Shears</p> <p>Designed for the rapid ligation and transection of tissue.</p>

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Accessory	Steris®	Sterrad®	EtO	Autoclave (Pre-vacuum)	Enzymatic Cleaner	Ultrasonic Bath
Type 1 Endoscopes (Straight and angled)		✓	✓		✓	
Type 2 Endoscopes (Straight and angled)	✓	✓	✓		✓	
Type 3 Endoscopes (Straight and angled)	✓	✓	✓		✓	
Light Guide Cable			✓ ^a	✓	✓	
Sterile Camera Arm Adapter				✓	✓	✓
Sterile Camera Adapter				✓	✓	✓
Sterile Instrument Arm Adapter				✓	✓	✓
Scope Alignment Targets				✓	✓	✓
Light Guide Cable Adapter				✓	✓	✓
Camera Cannula Mount				✓	✓	✓
Instrument Arm Cannula Mount				✓	✓	✓
Instrument Arm Cannula				✓	✓	✓
Endowrist Instruments				✓	✓	✓

a. The Bifurcated Light Guide is not compatible with EtO sterilization.

Process	Restricting Condition
EtO sterilization	Temperature: 55 ± 2°C Relative Humidity: 70 ± 5% Pressure Set Point: 25.4 PSIA Ethylene Oxide Concentration: 600 ± 30 mg/L Gas Exposure Time: 2 hours Detoxification Time: 0 hours Drying time: 0 hours Aeration: 12 hours at 55± 2°C

Autoclave Note: Do not autoclave endoscopes.	Cycle: Pre-vacuum Temperature: 270° - 272°F or 132° - 134°C Minimum Exposure Time: 4 minutes* Average Dry Time: 20 minutes
Sterrad® 100 Sterilization System	Vacuum Phase: 5-20 minutes Injection Phase: 6 minutes Diffusion Phase: 44 minutes Plasma Phase: 15 minutes Vent Phase: 5 minutes Total Cycle Time: Approximately 75 minutes
Steris System 1® Sterile Processing System	Sterilant: Peracetic Acid Sterilization time: 12 minutes Process temperature: 50 - 56°C Total processing time: 30 minutes Note: Please consult the Steris System 1® Sterile Processing System Operators Manual for detailed instructions for use.

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Internal System fault (reference: xxxx).

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Difficulty determining Arm position.

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1. Difficulty determining Arm position.
2. Arm movement not as expected.
3. Temporary motion sensor data mismatch.
4. Excessive force detected on Arm(s).
5. Power fluctuation detected.
6. System functioning with battery backup low.
7. Internal clock battery low.
8. System functioning without reserve power supply.

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符號	意義	位置	附註
一般使用			
	Attention: Consult accompanying documents	Several places on the da Vinci System	
	CF Applied Part	On instruments	Indicates instrument meets the requirements for a CF applied part.

符號	意義	位置	附註
一般使用			
	Equipotential Terminal	Rear panel of Surgeon Console, base of Surgical Cart and base of Vision Cart	Terminal for connecting Console or Cart to an equipotential point.
	Focus Controller	Focus Controller Box	Indicates where Focus Control Cable is inserted into the Focus Controller.
—	Focus Out	Surgeon Console Foot Switch Assembly and on Vision Cart Focus Controller	Moves focus of endoscope in negative direction.
+	Focus In	Surgeon Console Foot Switch Assembly and on Vision Cart Focus Controller	Moves focus of endoscope in positive direction.
	Flush Port	On instruments adjacent to Flush Port	Indicates connection point for flushing instrument.
S/N	Serial Number		


符號	意義	位置	附註
醫師控制中心			
	Composite Output	Rear panel of Surgeon Console	Low level composite video output (for future use).
	Focus	Rear panel of Surgeon Console adjacent to Focus Controller connector	Denotes where the Input/Output Focus Controller serial connector from the video cable is connected to the Console.
	Foot Switch	Rear panel of Surgeon Console	Indicates connector controls for Foot Switch function.
	Input and Output	Rear panel of Surgeon Console adjacent to Focus Controller connector	Denotes where the Input/Output Focus Controller serial connector from the video cable is connected to the Console.
	OFF for part of the da Vinci System	Rear panel of Surgeon Console adjacent to main circuit breakers	Circuit breaker indication.
	ON for part of the da Vinci System	Rear panel of Surgeon Console adjacent to main circuit breakers	Circuit breaker indication.
	Serial Digital Input	Rear panel of Surgeon Console	Low level, digital video input.
	Serial Digital Output	Rear panel of Surgeon Console	Low level, digital video output (for future use).
	Signal Input	Rear panel of Surgeon Console adjacent to ESU Foot Switch connector	Indicates input from Foot Switch of ESU.

符號	意義	位置	附註
醫師控制中心			
	S-Video Input	Rear panel of Surgeon Console	Low level S-video input (for future use).
	S-Video Output	Rear panel of Surgeon Console	Low level S-video signal output.
C	Chrominance	Rear panel of Surgeon Console	Low level video signal.
CAMERA	Camera Control	Surgeon Console Foot Switch Assembly	If the pedal is depressed and held down, the Masters disconnect from Instrument Arm control and connect to the Camera Arm.
CLUTCH	Clutch/Instrument Arm Swap pedal	Surgeon Console Foot Switch Assembly	If the pedal is depressed and held down, both Masters disconnect from all Instrument Arms.
COAG	Coagulation Function	Surgeon Console Foot Switch Assembly	Operable only if an appropriate compatible ESU is connected to the Console
EPO	Emergency Power Off	Surgeon Console Power Supply	Instantly disables the da Vinci System, including the control electronics and Stereo Viewer. CAUTION: To Power On after an EPO, turn off the main circuit breaker for 30 seconds.
L	Left	Rear panel of Surgeon Console	Indicates left Stereo Viewer signal for 3D video.

符號	意義	位置	附註
醫師控制中心			
R	Right	Rear panel of Surgeon Console	Indicates right Stereo Viewer signal for 3D video.
S1	Surgeon Console to Surgical Cart Cable 1	Rear panel of Surgeon Console	Indicates connection location for cable.
S2	Surgeon Console to Surgical Cart Cable 2	Rear panel of Surgeon Console	Indicates connection location for cable.
S3	Surgeon Console to Surgical Cart Cable 3	Rear panel of Surgeon Console	Indicates connection location for cable.
S4	Surgical Console to Surgical Cart Cable 4	Rear panel of Surgeon Console	Indicates connection location for cable.
SYSTEM	Instrument control system circuit breaker	Rear panel of Surgeon Console next to the da Vinci System circuit breaker	Controls main power to the da Vinci System only (not the Stereo Viewer).
VIEWER	Stereo Viewer circuit breaker	Rear panel of Surgeon Console next to Stereo Viewer circuit breaker	Controls main power to the Stereo Viewer only.
Y	Luminance	Rear panel of Surgeon Console	Low level video signal.

符號	意義	位置	附註
SURGICAL CART			
P1	Surgeon Console to Surgical Cart Cable 1	Rear of Surgical Cart	Indicates connection location for cable.
P2	Surgeon Console to Surgical Cart Cable 2	Rear of Surgical Cart	Indicates connection location for cable.
P3	Surgeon Console to Surgical Cart Cable 3	Rear of Surgical Cart	Indicates connection location for cable.
P4	Surgeon Console to Surgical Cart Cable 4	Rear of Surgical Cart	Indicates connection location for cable.
RIA	Connector designation for Intuitive Surgical Service personnel only	Non user-serviceable cables on the Surgical Cart	For Intuitive Surgical Service personnel only.

符號	意義	位置	附註
顯示器符號 (圖示)			
	Camera Control in progress (at Surgeon Console)	Stereo Viewer and Assistant Monitor	Appears when Camera/Scope is being positioned.
	Instrument Arm or Camera Arm clutching in progress	Stereo Viewer and Assistant Monitor	Appears while Instrument Arm or Camera Arm is being clutched.
	Instrument is invalid	Stereo Viewer and Assistant Monitor	Appears when Sterile Adapter is present and installed instrument is present and invalid.
	No uses left on Instrument	Stereo Viewer and Assistant Monitor	Instrument can no longer be used.
(flashing) 	Scope selection required at Surgeon Console	Stereo Viewer and Assistant Monitor	Camera/Scope combination must be selected after Homing.
(flashing) 	Instrument expires after procedure is complete	Stereo Viewer and Assistant Monitor	Instrument can be used during current procedure, but not for a new procedure.
(flashing) 	Instrument is being engaged	Stereo Viewer and Assistant Monitor	Amount of remaining use is displayed. In this example, nine uses remain for the installed instrument.

符號	意義	位置	附註
顯示器符號 (圖示)			
	No Sterile Adapter is Present	Stereo Viewer and Assistant Monitor	May also appear if Sterile Adapter is present, but not properly sensed by the da Vinci System.
	No instrument present	Stereo Viewer and Assistant Monitor	Appears while Sterile Adapter is present.
	Master repositioning in progress	Stereo Viewer and Assistant Monitor	Appears when Masters are being repositioned.
	Set-up Joint positioning	Stereo Viewer and Assistant Monitor	Appears while Set-up Joint is being positioned.
	Auxiliary Arm Z-axis movement	Stereo Viewer and Assistant Monitor	Appears while the auxiliary arm moves up or down along the Zaxis.
	Auxiliary Arm Z-axis Error	Stereo Viewer and Assistant Monitor	Appears when the movement of the Z-axis results in a safety check.
	The da Vinci System is operating on back-up back-up battery power. Battery is almost empty.	Stereo Viewer and Assistant Monitor CAUTION: During a general system power failure, the Stereo Viewer does not receive power. The icon will appear on the Assistant Monitor only.	Appears when there is a general power failure or when there is an internal system power supply failure and the battery has been run long enough to consume almost all of its energy.

	<p>The da Vinci System is operating on back-up battery power. Battery is full.</p>	<p>Stereo Viewer and Assistant Monitor</p> <p>CAUTION: During a general system power failure, the Stereo Viewer does not receive power. The icon will appear on the Assistant Monitor only.</p>	<p>Appears when there is a general power failure or when there is an internal system power supply failure and the battery is fully charged.</p>
	<p>The da Vinci System is operating on battery power. Battery is half full.</p>	<p>Stereo Viewer and Assistant Monitor</p> <p>CAUTION: During a general system power failure, the Stereo Viewer does not receive power. The icon will appear on the Assistant Monitor only.</p>	<p>Appears when there is a general power failure or when there is an internal system power supply failure and the battery contains half of its available energy.</p>
	<p>Match Grip</p>	<p>Stereo Viewer and Assistant Monitor</p>	<p>Appears at the top left or right of the screen, next to the associated MTM. User must match the grip of the Master to the corresponding. Instrument tip to enter Follow Mode.</p>
	<p>Homing Failed</p>	<p>Stereo Viewer and Assistant Monitor</p>	<p>Appears when the homing procedure has failed.</p>

PSMs, ECMs and MTMs each have a portion of the screen dedicated to their icons. The Match Grip and Homing Failed icons appear in the dedicated screen area for each manipulator.

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Chapter 9 Glossary of Terms

2D	Two-dimension.
3D	Three-dimension.
3D Display	The three dimensional image created by two cameras. The image is viewed by looking into the Stereo Viewer on the Surgeon Console.
3D Viewer	See Stereo Viewer.
ACU	Auxiliary Control Unit. Adapter See Sterile Adapter.
Anchor Pads	Refers to the brakes on the Surgical Cart.
Arms	The part of the da Vinci System that holds a grip, instrument or camera. There are up to five arms on the da Vinci System. Two Master arms with Master grips are located on the Surgeon Console.
Console.	The Surgical Cart has up to three arms; one for the camera and two for the instruments.
Assistant Monitor	The monitor on the Vision Cart that displays a 2D image.
Auxiliary Instrument Arm	Some systems are available with a third instrument arm, also referred to as the "fourth arm." The auxiliary arm can be assigned to the left or right Master (or no Master). The arm can be activated with a pedal-click on the Clutch foot switch.
Camera Arm	The arm on the center Set-up Joint of the Surgical Cart that controls the camera/endoscope according to the surgeon's movements of the Masters.
Camera Orientation Line	Shows the orientation of the camera, with respect to the floor, in the lower right-hand corner of the Stereo Viewer and Assistant Monitor.

Cannula Mount	The accessory that attaches a cannula to an arm. There are two types: instrument arm cannulae and camera arm cannulae. Two styles exist for the instrument arms (standard and Twist-lock) and one style for the camera arm.
Cart	Refers to the Surgical Cart.
Carriage	The portion of the instrument arm to which an Instrument Sterile Adapter attaches. The carriage moves up and down the arm under control of the system, or manually by using the instrument arm clutch.
CRT	Cathode Ray Tube.
CAUTION	An important level of concern. Failure to heed a CAUTION may result in unintended motion of the da Vinci System that may result in injury to a patient.
Clutch (verb)	<p>A. Master clutch is the act of disengaging the Masters from the instrument arms and camera arm so the Masters can be repositioned in a more comfortable working space for the surgeon. This action is similar to lifting a computer mouse off the mouse pad and repositioning it. A clutch is also used at the Surgical Cart to position the instruments and the camera/endoscope.</p> <p>B. A Patient-Side clutch allows the patient-side manipulator to float, allowing manual adjustment of the instrument or camera arms.</p>
CLUTCH (noun)	<p>A. Foot Switch pedal used to control Master Clutching, and as part of the pedal-click to swap control of an instrument arm to the auxiliary instrument arm, if available.</p> <p>B. PSM or ECM button that allows clutching of the PSM or ECM.</p>
COAG	Foot Switch pedal to control Coagulation.
Console	See Surgeon Console.
DANGER	The highest level of concern. Failure to heed a DANGER warning can result in injury to a patient.
Drape Clip	A reusable accessory that reinforces the interface between the adhesive strip on the instrument arm drape and the instrument arm Sterile Adapter.

ECM	Endoscopic Camera Manipulator. This acronym refers to the camera arm.
ESU	Electrosurgical Unit or Electrosurgical Generator Unit.
Endoscope	An instrument used for the examination of the interior of a canal or a hollow space; also called a "scope."
Endowrist® Instruments	Instruments with a wrist located near the tip.
Emergency Stop State	A state where the motors of are turned off and a screen message is sent to the user.
Faulting	The transition from a working state to either a soft-locked or dampened state.
Foot Switch	A pedal or switch located on the Foot Switch Assembly that has a specific function.
Foot Switch Assembly	Located at the base of the Surgeon Console, containing all of the foot controls.
Head Sensor	Infrared sensor on either side of the View Port of the Surgeon Console located above the Stereo Viewer.
Illumination System	See Light Source.
Instrument	Any one of several tools used to effect the procedure in the patient once attached to an instrument arm and inserted into the patient. Instruments include, for example, Large Needle Drivers, DeBakey Forceps and Round Tip Scissors.
Instrument Arms	The arms on the outer Set-up Joints of the Surgical Cart that manipulate the instruments according to surgeon's movement of the Masters.
LED	Light emitting diode.
Light Source	Endoscopic illumination system. In the da Vinci System, there are two light sources. Together they provide illumination inside the body for vision.
Master	The control arms and grips in the Surgeon Console that the surgeon grasps and moves. The surgeon's movements are translated to the instruments and camera attached to the Surgical Cart arms.
MIS	Minimally Invasive Surgery.

MTM	Master Tool Manipulator. Refers the Masters on the Surgeon's Console.
Notes	User information to emphasize an important point.
Pedal Click	The act of quickly tapping a foot switch at the Surgeon's Console to alter configuration. The pedal click is used with the Clutch foot switch to swap control between instrument arms, and with the Camera foot switch to toggle between Narrow (3D) and Wide (2D) angle views.
PSM	Patient-side manipulator (instrument arm and respective Set-up Joint).
READY	A button on the User Switch Panel of the da Vinci System used to change the da Vinci System state.
Scope	Endoscope.
Screen	The monitor display, located on the Stereo Viewer and/or the Assistant Monitor.
Setup Joint	The joints on the Surgical Cart that support the instrument arms and the camera arm. These joints are used to set up the initial positions of the arms on the Surgical Cart.
STANDBY	A button on the USP to put the da Vinci System into the Stand By mode.
Stereo Viewer	The fixed eyepiece located on the Surgeon Console where the surgeon looks to see the 3D image.
Sterile Adapter	Interface device that allows the sterile barrier to be maintained between the camera arm and the endoscope or the instrument arm and the instrument attached to the arm. There are various types of Sterile Adapters: a Sterile Camera Adapter, Sterile Camera Arm Adapter and an Instrument Sterile Adapter. They are not interchangeable.
Sterile Adapter Reinforcement	An adhesive strip on the drape that is adhered to the Sterile Adapter.
Surgeon Console	The part of the da Vinci System consisting of a structure that supports the Masters and the 3D Viewer. The Surgeon Console also houses most of the electronics of the complete da Vinci System and the power supply and battery.

Surgical Cart	The part of the da Vinci System that is located on the patient side and consists of the column that supports the Set-up Joints that in turn support the instrument and camera arms.
Surgical Cart Arms	The arms that are components of the Surgical Cart: two or three instrument arms and one camera arm, depending on the system.
UIP	User Interface Panel. The set of buttons and LEDs located on the left side of the Surgeon's Console, as you face it.
User Interface	The controls located on either side of the armrest of the Surgeon Console that are used to command the da Vinci System. The User Interface consists of the User Interface Panel and the User Switch Panel.
USP	User Switch Panel. The set of buttons located on the right side of the Surgeon's Console, as you face it.
View Port	The area on the Surgeon Console that includes the View Port, Stereo Viewer and infrared Head Sensors.
WARNING	High level of concern. Failure to heed a warning could result in harm to a patient.

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Accessories

Cleaning 7-1

Sterilization and Cleaning matrix 7-2

ALARM indicator 6-5

Alarms

Audible alarms 8-3

Alignment Target

Using for Calibration 3-9

Audio Intercom System

Volume adjustment 6-11

Autoclave

Sterilization reference information 7-3

Auxiliary Instrument Arm

Assigning to a master 4-4

Collision Detection

Reference information 6-18

Disabling 6-8

Display information

Stereo Viewer

Auxiliary Instrument Arm information

6-8

Excessive Force Detection

Reference information 6-19

Moving along the Z-axis 6-17

Reference information 6-16

Swapping 4-5

Reference information 6-14

Using 4-4

B

Blade Protector 4-12

C

Calibrating the Endoscope and Camera 3-9
Calibration
TEST IMAGE button
Reference information 6-7
Using the TEST IMAGE button 3-9
Camera
Adjusting with the Camera Orientation Line
4-4
Attaching to the Camera Arm 4-17
Calibrating 3-9
Repositioning the Camera/Scope image 4-3
Camera Arm
Attaching the Camera Cannula Mount 3-6
Attaching the Sterile Camera Arm Adapter
3-7
Attaching/Removing the Camera/Endoscope
4-17
Cannula connection 4-7
Clutching 4-8
Draping 3-5
overview 1-5
Reference information 6-16
Removing the Cannula Mount 5-1
Removing the Sterile Camera Adapter 5-1
Camera Cannula Mount
Removing from the Camera Arm 5-1
CAMERA Foot Switch
Using to reposition the Camera/Scope image
4-3
Camera Head 3-7
Attaching to the endoscope 3-8
Camera Orientation Line
Using to adjust the camera 4-4

Camera Signal Loss indicator 6-4
Camera Sterile Adapter
Attaching to the endoscope 3-7
CANNULA button
Reference information 6-6
Selecting the Cannula Length 4-6
Cannula Mount
Attaching the Camera Cannula Mount 3-6
Attaching to instrument arms 3-4
Removing from the Instrument Arms 5-1
Cannulae
Attaching the Camera Cannula Mount 3-6
Attaching the Cannula Mount 3-4
Attaching to Instrument arms 3-4
Attaching Twist-lock Cannulae 3-4
Connecting the Camera Cannula to the Camera
Arm 4-7
Connecting to the Instrument Arms 4-7
Inserting and Connecting 4-6
Selecting the Cannula Length 4-6
Cleaning 7-1
Accessories 7-1
Sterilization and Cleaning matrix 7-2
Surgeon's Console 7-1
Surgical Cart 7-1
Clip Appliers
Index-2 Intuitive Surgical, Inc.
P/N 550120-02 Rev XB
Loading and firing 4-15
CLUTCH Foot Switch
Using to reposition Masters 4-5
Using to swap instrument arms 4-5
Clutching

Camera clutching 6-15
Instrument and Camera arm clutching 4-8
Master Clutching
Reference information 6-14
CRITICAL ERROR indicator 6-5

D

da Vinci System
Draping 3-1
General Overview 1-1-1-2
General setup procedures 2-1-2-6
Shutdown 5-1
Startup 2-5-2-6
Use during a procedure 4-1-4-17
DISPLAY button
Reference information 6-6
Draping 3-7
Attaching sterile adapters 3-2
Attaching the Sterile Camera Adapter 3-7
Camera Arm
3-5
Draping the Camera Head 3-7
Instrument arms 3-2
Reference information 6-19
Securing the drapes 3-5
Sequence 3-2

E

Emergency Grip Release 4-13
Emergency Power Off button
Reference information 6-11
Emergency Stop
reference information 8-1
EMERGENCY STOP button 6-3
Emergency Stop mode 3-12

Endoscope

Attaching the Sterile Camera Adapter 3-7

Attaching to the Camera Arm 4-17

Attaching to the Camera Head 3-8

Calibrating 3-9

Repositioning the Camera/Scope image 4-3

Selecting the scope angle 4-2

Using the Alignment Target 3-9

Error Log

Viewing 8-4

Errors

Fault and Emergency Stop states 8-1

EtO

Sterilization reference information 7-2

F

FAULT OVERRIDE button 6-5

Faults

Loss of Power 8-3

Non-recoverable 8-2

Recoverable 8-2

reference information 8-1

Focusing the video image 4-3

Follow mode 3-11

Foot Switches

Camera 6-14

Clutch 6-14

Focus 6-15

G

Guided Tool Change 4-10

H

Head Sensor

Reference information 6-10

Homing 2-5

I

Icons

Reference information 6-12, 8-4–8-11

InSite Vision System

General overview

Vision System

General Overview 1-6

Instrument Arms

attaching cannulae 3-4

Clutching 4-8

Connecting cannulae 4-7

Draping 3-2

overview 1-4

Reference information 6-16

Removing sterile adapters 3-3

Removing the Cannula Mount 5-1

Intuitive Surgical, Inc. Index-3

P/N 550120-02 Rev XB

Removing the Sterile Adapters 5-1

Sterile Adapter compatibility 3-13

Swapping 4-5

Instruments 7-1

Automated insertion 4-10

Cleaning 7-1

Clip Appliers

Loading and firing 4-15

Emergency grip release 4-13

Limited Reuse and Expiration indication

6-22

Manual insertion 4-10

Placement and removal 4-9

Reference information 6-19–6-22

Removal 4-11

Scalpel blade insertion instructions 4-14

Tool Reader button 6-10

Tool Reader mode 4-14

Usage summary display 5-2

Using a blade protector 4-12

L

Loss of Power 8-3

M

Main circuit breakers 2-2

Masters

Assigning to Auxiliary Instrument Arm 4-4

Assigning to Instrument Arms

Reference information 6-7

Clutching 6-14

General overview 1-2

Repositioning during use 4-5

Swapping 4-5

Using to reposition the Camera/Scope image

4-3

Modes

Emergency Stop 3-12

Follow mode 3-11

Ready mode 3-10

Stand By mode 3-10

Understanding Operating Modes 3-10–3-12

O

Operating Modes 3-10–3-12

P

Port Locations

Selecting 4-1

Power Loss indicator 6-4

R

READY button 2-6

Reference information 6-4

Ready mode 3-10

Remote Center

Reference information 6-16

S

SCALING button

Using to change scaling options 4-5

SCALING control

Reference information 6-5

Scalpel blades

Insertion instructions 4-14

SCOPE ANGLE button

Reference information 6-5

using for calibration 3-9

SELECT MASTER buttons

Reference information 6-7

Select Master buttons

Using to assign a master to an instrument
arm 4-4

Selftest 2-5

Service Required indicator 6-4

Setting up the da Vinci System 2-1–2-6

Setup Joint

overview 1-4

Reference information 6-15

Shutdown process 5-3

Shutting down the system 5-1

SILENCE ALARM button 6-5

STAND BY button

Reference information 6-4

Stand By mode 3-10

Starting up the system 2-5–2-6

Startup

Homing the arms 2-5

READY button 2-6

Selftest 2-5

Start-up Sequence indicator 6-4

Stereo View

Index-4 Intuitive Surgical, Inc.

P/N 550120-02 Rev XB

Brightness and Contrast adjustment 6-11

Stereo Viewer

DISPLAY button 6-6

Focusing the video image 4-3

General Overview 1-2

Message and Icon locations

Icons

Display locations 8-12

Switching from 3D to 2D 4-3

Sterile Accessories

Example 3-1

Preparing 3-1

Sterile Adapter compatibility matrix 3-13

Sterile Adapters

Attaching to drapes 3-2

Reference information 6-19

Removing 5-1

Removing from Instrument arm 3-3

Sterile Camera Adapter

Attaching to the Camera Head 3-8

Removing from Camera Arm 5-1

Sterile Camera Arm Adapter

Attaching to Camera Arm 3-7

Sterilization

Autoclave

Reference information 7-3
Steris
Reference information 7-3
Sterrad
Reference information 7-3
Sterilization
EtO
Reference information 7-2
Steris
Sterilization reference information 7-3
Sterrad
Sterilization reference information 7-3
Surgeon's Console 7-1
Adjustment for comfort 4-2
Audio Intercom System volume adjustment
6-11
Brightness and Contrast adjustment 6-11
Controls and Indicators
General overview 1-2–1-3
Reference information 6-1–6-15
Emergency Power Off button 6-11
Focusing the video image 4-3
Foot Switches 6-14
General Overview 1-2
Head Sensor
Reference information 6-10
Intra-ocular Spacing Adjustment 6-11
Masters
General overview 1-2
Repositioning the Camera/Scope image 4-3
Selecting the endoscope angle 4-2
Setup procedures 2-1–2-2
Main breakers 2-2

Positioning 2-1
Startup 2-5–2-6
Stereo Viewer
General overview 1-2
Switching from 3D to 2D 4-3
User Interface Panel 6-2
User Switch Panel 6-2
Viewer Height adjustment 6-10
Surgical Cart
Camera Arm
overview 1-5
Reference information 6-16
Cannula Insertion and Connection 4-6
Cleaning 7-1
General overview 1-4–1-5
Instrument Arms
overview 1-4
Reference information 6-16
Setup
Connecting system cables 2-2
Positioning 2-2
Setup Joint
overview 1-4
Reference information 6-15
Startup 2-5–2-6
Symbols
Reference information 8-4–8-11
SYSTEM button
Reference information 6-3
Using for shutdown 5-3
System Status indicators
Alarm 6-5
Critical Error 6-5

Fault Override 6-5

Silence Alarm 6-5

Intuitive Surgical, Inc. Index-5

P/N 550120-02 Rev XB

T

TEST IMAGE button

Reference information 6-7

Using for calibration 3-9

TOOL READER button

Reference information 6-10

Tool Reader mode 4-14

Twist-lock Cannulae

Attaching to instrument arms 3-4

U

Usage Summary Display 5-2

User Interface Panel

Example 6-2

READY button 6-4

STAND BY button 6-4

Status Indicators

Camera Signal Loss 6-4

Power Loss 6-4

Service Required 6-4

Start-up Sequence 6-4

User Switch Panel

EMERGENCY STOP button 6-3

Example 6-2

SYSTEM button 6-3

V

Video image

repositioning 4-3

Vision Cart

Setup

Connecting the Vision Cable 2-3

Positioning 2-3

W

Wide Angle View

Toggling between 3D and 2D 6-15

Z

Z-axis

Moving the Auxiliary arm 6-17

Index-6 Intuitive Surgical, Inc.

P/N 550120-02 Rev XB